

▼
**NECK
MASSAGER**
▲



www.belmint.com



support@belmint.com



1 (844) 423-5646



@belmintproducts



@belmintproducts

**NECK
MASSAGER**



▶ User Guide

 Model	 Voltage	 Power Consumption
BEL-BNM	12V	24W

Congratulations on your purchase!

Your tootsies are in for a treat!

Your Belmint Neck Massager shiatsu neck & back massager with heat makes it possible for you to get a spa-level massage anywhere, anytime. This compact unit comes in its own custom carry case, which makes it easy to transport. Now you can experience professional spa-level massage at work, in the car (not while driving please!) or on a plane.

Like all Belmint products your portable neck & back massager is crafted with high-quality, durable materials to ensure your comfort. Specially designed to ease tension in tired and sore muscles, your massager is equipped with rolling and kneading nodes that enables you to target those hard-to-reach areas of your body.

For a truly custom massage experience your portable neck & back massager is equipped with optional heat, directional control and nodular speed control—because some days we need more kneading than others. Cleanup and storage are easy. Simply wipe with a damp cloth after use and store it in its included carry case. You won't want another portable neck & back massager. You've got a Belmint!

To celebrating your comfort,

The Belmint Team

▶ Sometimes peace is hard to find

But good service shouldn't be

At Belmint, we strive to provide our customers with ultimate relaxation.

If you encounter any issues with the performance of your device, our authorized service personnel are here to help. Check back cover for information.

▶ What's in the Box?

- 1 Neck massager with Heat
- 1 DC Plug
- 1 AC Plug
- 1 User Manual Guide

▶ How to Use

1. The power adapter should be connected with the power socket. If you're using the massager in the car, plugged in the power cord to the car cigarette lighter socket.
2. Press the button on the controller to turn on the massager.
3. Use the controller option for your desired massaging options.
4. After use, press the power button in the controller to turn off the massager and unplug the power adapter from the socket or power supply.